

410-VIS-IR REFLECTOMETER

Solar reflectance and thermal emittance measurements

The 410-Vis-IR joins the ET100 and 410-Solar or 410 Solar-i measurement heads in a single package for measuring solar reflectance/absorptance and thermal emittance. The ET100 measurement head collects directional reflectance from 1.5 to 21 μm to determine directional and total hemispherical emissivity. The 410-Solar models measure reflectance at near-normal incidence from 335-2500 nm. Easily switch between reflectance and emittance data collection with a single handle for data display and measurement controls.



BENEFITS

- Complete spectral range**
 Reflectance, absorptance and emittance from 350 nm – 21 μm .
- ASTM compliant**
 Compliant with ASTM E903, E408, C1549 & E1980.
- No waiting**
 90 sec warm up, no equilibration between measurements.
- Easy to use**
 Simply press port against sample and pull trigger.
- Fast calibration**
 One minute calibration at start of measurement session.
- Room temperature samples**
 Calculate emissivity without heating sample.

APPLICATIONS

- Space Coatings**
 Thermal control | α/ϵ | Thermo-optical properties
- Defense | Aerospace**
 IR Signature | Low observable paint & coatings
- Cool Building Materials**
 TSR | SRI | ASTM | LEED | CRRC
- Concentrated Solar**
 Mirror evaluation | Selective absorber coatings
- Radiative Heat Transfer**
 Absorptance for thermal modeling
- Semiconductors**
 Wafer fab hardware emissivity
- Astronomy**
 Mirror evaluation

ORDERING

Standard components	0410-0011	410-Vis-IR Measurement Heads
410-Solar	0410-0001	Handheld Command Module - 120VAC
	0410-0100	Specular Gold Calibration Coupon (Non-NIST Traceable)
	0410-0102	Diffuse Calibration Coupon (Non-NIST Traceable)
	0410-0103	Specular Calibration Coupon (Non-NIST Traceable)
Standard components	0410-0013	410-Vis-IR (Model i) Measurement Heads
410-Solar-i	0410-0001	Handheld Command Module - 120VAC
	0410-0100	Specular Gold Calibration Coupon (Non-NIST Traceable)
	0410-0107	Glazed Ceramic Calibration Coupon (Non-NIST Traceable)
Options	0410-0002	Benchtop Remote Control Unit - 120VAC
	0410-0101	Specular Gold Calibration Coupon (NIST Traceable)
	0410-0106	Diffuse Calibration Coupon (NIST Traceable)
	0410-0104	Specular Calibration Coupon (NIST Traceable)
	0410-0123	Glazed Ceramic Calibration Coupon (NIST Traceable)
	0410-1015	410-Vis-IR Reflectometer Maintenance and Calibration Plan (Non-NIST)
	0410-1014	410-Vis-IR Reflectometer Maintenance and Calibration Plan (NIST)
	0410-1007	410-Vis-IR Extended Warranty
	0410-1012	410-Vis-IR (Model i) Extended Warranty
	0410-0200	Handheld Command Module - 220VAC
	0410-0019	Benchtop Remote Control Unit - 220VAC

EXAMPLE MENU SCREENS

Measurement screen. Results are displayed on the liquid crystal display touchscreen, and stored on a SecureDigital (SD) card.

nm	Specular	Total
335-380	.925±.001	.926±.001
400-540	.940±.002	.942±.002
480-600	.947±.001	.950±.001
590-720	.948±.001	.953±.001
700-1100	.950±.002	.955±.002
1000-1700	.951±.001	.956±.002
1700-2500	.953±.001	.958±.001

Solar absorptance calculation for the selected solar irradiance function.

Solar Irradiance:
Air Mass 0

Solar Absorptance:
 $\alpha_t = .116 \pm .001$
 $\alpha_s = .114 \pm .001$

Directional and hemispherical emittance measurement data screen.

Thermal Emittance

F:1 S:1
T:300 K Dielectrics

$\epsilon_{20^\circ} = 0.612 \pm 0.001$
 $\epsilon_{60^\circ} = 0.778 \pm 0.001$
 $\epsilon_H = 0.604 \pm 0.001$

Set temperature and material type for Hemispherical Thermal Emittance calculation.

Thermal Emittance

engmeas1/sample1
Temperature [K]: 300.0
HTE table: dielectrics

DTE at 20 deg: 0.612
DTE at 60 deg: 0.778
HTE: 0.604

CE Equipment described herein is subject to US export regulations and may require a license prior to export. Diversion contrary to US law is prohibited.

SPECIFICATIONS 410-SOLAR, 410-SOLAR-i

410 Solar or 410 Solar-i included in 410 Vis-IR Package.


	410-SOLAR	410-SOLAR-i
MEASURED DATA		
<i>Measured Parameter</i>	Directional hemispherical reflectance (DHR)	
<i>Method</i>	Integrated total reflectance in a band for a given angle of incidence	
<i>Measured Value</i>	Absolute reflectance (0-1), Diffuse Reflectance	Absolute reflectance (0-1)
<i>Calculated Value</i>	Total Solar reflectance, Solar absorptance, specular reflectance	Total Solar reflectance, Solar absorptance
<i>Wavelength Bands (nm)</i>	335-380, 400-540, 480-600, 590-720, 700-1100, 1000-1700, 1700-2500	
<i>ASTM Standards</i>	C1549 E1980 E903	
<i>Angle of Incidence</i>	20° from normal incidence	
<i>Calibration Coupon</i>	Specular, Diffuse	Glazed Ceramic
PERFORMANCE		
<i>Accuracy</i>	+/- .02	
<i>Repeatability</i>	±.005 units	
<i>Beam Spot Size</i>	0.25 inches	0.50 inches
<i>Beam Angle</i>	3° half cone angle	N/A
<i>Measurement Time</i>	10 seconds	7 seconds
<i>Solar Irradiance Functions</i>	Air Mass 0 (AM) Extraterrestrial irradiance (ASTM E490-00) Hazy sky AM1.5 beam-normal irradiance (ASTM E891-87) Clear sky AM1 global horizontal irradiance (SMARTS 2.9.5) Clear sky AM1.5 global irradiance surface tilted 37° (ASTM G173-03) Clear sky AM1.5 global irradiance surface tilted 20° (ASTM G197-14) Clear sky AM1.5 global irradiance surface tilted 90° (ASTM G197-14) Clear sky AM1.5 global horizontal irradiance (SMARTS 2.9.5) Clear sky AM2.0 global horizontal irradiance (SMARTS 2.9.5)	
<i>Sample Size and Geometry</i>	Flat: ≥ 0.5 in. diameter Curved: 6 in. convex; 12 in. concave	
<i>Warm Up Time</i>	90 seconds	
<i>Time Between Measurements</i>	2 seconds	
<i>Sample Temperature</i>	Ambient or heated/cooled to 0 - 100° C	
<i>Operating Temperature</i>	0° to 40° C	
POWER		
<i>Run Time</i>	2 hours on one battery. Battery easily replaced with continuous operation after battery replacement.	
<i>Power Source</i>	Rechargeable battery (standard environmentally friendly NiMH)	
<i>Battery Recharge Time</i>	1 hour	
<i>VIS-NIR Source</i>	Tungsten filament, temperature controlled by user	
DIMENSIONS		
<i>Weight</i>	4.7 lbs. (2.13 kg)	
<i>Form Factor/Size</i>	H 11.54 in., L 9.04 in., W 3.27 in. (29.31 cm x 22.96 cm x 9.44 cm)	
INTERFACE		
<i>Operator Interface</i>	LCD graphics screen, 1/4 VGA, touch screen, software buttons; trigger switch in handle	
<i>Diagnostics</i>	On screen status and signals monitor. Signal values stored with data. Raw data collection and display.	
MISCELLANEOUS		
<i>Format</i>	Data files can be opened and post processed with Excel or a text processor	
<i>Storage</i>	Removable SanDisk (SD) card	
<i>Export control</i>	ECCN #3A999.F	

SPECIFICATIONS ET-100

ET-100 included in 410 Vis-IR Package.

ET-100	
MEASURED DATA	
<i>Measured Parameter</i>	Directional hemispherical reflectance (DHR)
<i>Method</i>	Integrated total reflectance in a band for a given angle of incidence
<i>Measured Value</i>	Absolute reflectance (0-1)
<i>Calculated Value</i>	Directional thermal emissivity at 20°, directional thermal emissivity at 60°, hemispherical thermal emissivity
<i>Wavelength Bands (microns)</i>	1.5-2.0, 2.0-3.5, 3.0-4.0, 4.0-5.0, 5.0-10.5, 10.5-21
<i>Angle of Incidence</i>	20° & 60° from normal incidence
<i>ASTM Standards</i>	E903
<i>Calibration Coupon</i>	Specular Gold
PERFORMANCE	
<i>Accuracy</i>	+/- .03
<i>Repeatability</i>	±.005 units
<i>Beam Spot Size</i>	0.50 inches
<i>Measurement Time</i>	10 sec
<i>Sample Size & Geometry</i>	Flat: ≥ 0.5 in. diameter Curved: 6 in. convex; 12 in. concave
<i>Warm Up Time</i>	90 seconds
<i>Time Between Measurements</i>	2 seconds
<i>Sample Temperature</i>	Ambient or heated/cooled to 0 - 100° C
<i>Operating Temp</i>	0° to 40° C
POWER	
<i>Run Time</i>	2 hours on one battery. Battery easily replaced with continuous operation after battery replacement.
<i>Power Source</i>	Rechargeable battery (standard environmentally friendly NIMH)
<i>Battery Recharge Time</i>	1 hour
<i>IR Source</i>	Kanthal filament operated at about 1,000°C
ENVIRONMENT	
<i>Storage</i>	-25° to 70°C
<i>Operating</i>	0° to 40° C
DIMENSIONS	
<i>Weight</i>	4.7 lbs
<i>Form Factor/Size</i>	H 11.54", L 9.04", W 3.27" (29.31 cm x 22.96 cm x 9.44 cm)
INTERFACE	
<i>Operator Interface</i>	LCD graphics screen, 1/4 VGA, touch screen, software buttons; trigger switch in handle
<i>Inspection Applications</i>	Pass/fail can be incorporated, user set values
<i>Diagnostics</i>	On screen status and signals monitor. Signal values stored with data. Raw data collection and display.
MISCELLANEOUS	
<i>Date Format</i>	Data files can be opened and post processed with Excel or a text processor
<i>Data Storage</i>	Removable SanDisk (SD) card
<i>Export control</i>	ECCN #3A999.F

410 SERIES REFLECTOMETERS & EMISSOMETERS




METHODOLOGY

The basic structure of a measurement head is an internal source, a modified integrating sphere, and detectors. The reflectance measurement is made by collimating the source beam onto the target, the energy is reflected back into the sphere, and eventually detected or dissipated.

The 410 Series Reflectometers measures the integrated surface reflectance of a surface at a given angle of incidence (20° or 60°). The integrating sphere captures the reflected light from the target material, integrating reflections in all directions. Wavelength-filtered detectors measure the total light reflected in each wavelength band and converts it to an analog electrical signal.

The 410 Series Reflectometer electronics processes the detector signals for initial amplification (fixed), filtering, offset adjustment, secondary amplification (variable), and analog to digital conversion. The digitized signals are read by the on-board processor, stored in memory, and then used to determine the target sample reflectance at each incident angle and wavelength band. Those reflectances are used to calculate additional properties such as directional thermal emittance or total hemispherical emittance. Results are displayed on the liquid crystal display touchscreen, and stored on a SecureDigital (SD) card.



INTEGRATING SPHERE SCHEMATIC

Schematic of the integrating sphere in contact with a sample.
Red arrow - illuminating beam
Purple arrow - reflected beam
Green arrows - scattered light

CALIBRATION COUPON



CALIBRATION

An easy calibration process is required before each measurement session. The software GUI will walk the user through the process. Calibration is performed using calibration coupon(s) with known reflectance values.

410 REFLECTOMETERS MODEL COMPARISON GUIDE

The SOC410 Series Reflectometers are portable contact measurement devices designed to take precise, accurate reflectance and emittance measurements. Made with an ergonomic power-drill design, the SOC410 Series lets you easily take measurements in-the-field or around the lab—no cords or external batteries necessary. The world's largest defense, aerospace, and energy companies rely on SOC410 data.



Model	410-Solar	410-Solar-i	410-VIS-IR	ET-100	ET-10	410-DHR
<i>Spectral Bands</i>	335 - 380 nm 400 - 540 nm 480 - 600 nm 590 - 720 nm 700 - 1100 nm 1000 - 1700 nm 1700 - 2500 nm	335 - 380 nm 400 - 540 nm 480 - 600 nm 590 - 720 nm 700 - 1100 nm 1000 - 1700 nm 1700 - 2500 nm	Dual measurement head package consisting of a 410-Solar model and ET100 measurement heads with a single command module	1.5 - 2.0 μm 2.0 - 3.5 μm 3.0 - 4.0 μm 4.0 - 5.0 μm 5.0 - 10.5 μm 10.5 - 21.0 μm	3.0-5.0 μm 8.0-12.0 μm	0.9 - 1.1 μm 1.9 - 2.4 μm 3.0 - 4.0 μm 3.0 - 5.0 μm 4.0 - 5.0 μm 8.0 - 12.0 μm
<i>Calculated Properties</i>	Total, diffuse & specular reflectance absorptance	Total reflectance/absorptance		In-band total reflectance Directional thermal emissivity at 20° Directional thermal emissivity at 60° Hemispherical thermal emissivity	Directional thermal emissivity at 20°	In-band total reflectance In-band emissivity
<i>Angle of Incidence</i>	20°	20°		20° and 60°	20°	20° and 60°
<i>Calibration Coupon(s)</i>	Solar Diffuse Solar Specular	Glazed Ceramic		Specular Gold	Specular Gold	Specular Gold
<i>ASTM Compliance</i>	C1549 E903 E1980	C1549 E903 E1980		E408 E1980		N/A



SURFACE OPTICS CORPORATION

11555 Rancho Bernardo Road
San Diego, CA 92127

For Information and Ordering

Email: contact@surfaceoptics.com
Phone: +1 858 675-7404
Website: surfaceoptics.com

This information is subject to change without notice.
© Surface Optics Corporation 2021